
Instructions

Storage up
to 30 °C

18 m

Sample Preparation

Viable count in solid foodstuff

Take 25 g of solid specimen into the sterilized homogenizer bag. Add 225 mL of sterile Buffered Peptone Water or EEM Broth into the bag, and homogenize with stomacher for about one (1) minute. Incubate the bag 20-24 hours at 35–37 °C to pre-enrich.

Viable count in water/liquid food

Add 9 part of buffered peptone water or EEM broth to 1 part of liquid sample. Filter the liquid sample through membrane filter, and put the filter into BPW or EEM Broth. Incubate the bag 20-24 hours at 35–37 °C in the incubator for pre-enrichment.

Viable count in swab test specimen

Add 1mL of wiping solution to 9 mL of peptone water or EEM broth.

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6

CompactDry Protocol

1 Open the cap **2** Dispense 0.1 mL of enriched sample on the plate (approx. 1 cm from the plate's rim). **3** Dispense 1 mL of sterilized water at the opposite point. The water will diffuse uniformly across the plate. **4** Replace the cap, label and turn the plate over. **5** Incubate for the required time and temperature. **6** Following incubation, count the number of coloured colonies.

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Incubation temperature: 41 -43 °C
Incubation time: 20 -24 hours

Interpretation

Salmonella Positive - Black to green isolated or fused colonies. The medium around the colonies changes to yellow. If a large number of Salmonella are present, the medium becomes completely yellow.

Salmonella Negative - No yellow colour within the medium. The medium becomes reddish purple colour and no black or green colonies observed.

General information

Caution: Plate may become yellow due to Pseudomonas or Proteus, but yellow portion is small due to lower motility.

Final report for Salmonella positive or negative result requires identification and confirmation tests. Salmonella can be isolated from both black/green colonies and/or from the yellow portion away from the inoculation point. Isolation of Salmonella from colonies: Take black to green colonies with loop and inoculate on an appropriate agar (e.g. MLCB). Post-culture, continue with conventional confirmation testing. After use please follow local disposal regulations.

Storage and shelf life

Store at room temperature (+1 to +30 °C). Shelf life 18 months after manufacturing.

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