

Accu-Bac

Quantitative Micro-Organisms



About

Accu-Bac is a quality control micro-organism. It is a guaranteed second generation derivative of NCTC[®] and NCPF[®] cultures, manufactured under licence from UK Health Security Agency. Accu-Bac is to be used in an industrial setting. Each carefully prepared disc contains a target range of CFU, with 95% confidence interval stated on the Certificate of Analysis. This is ideal for most quantitative applications. Accu-Bac is supplied in 1 disc vials and has at least a 12 month shelf-life from the date of manufacture.

This product is not currently CE marked.



TCS Biosciences Ltd
+44 (0)1296 714222
sales@tcsgroup.co.uk
www.tcsbiosciences.co.uk



Accu-Bac

Quantitative Micro-Organisms



- Defined target range CFU counts per disc, with 95% confidence interval stated on the Certificate of Analysis
- Easy to use by following the simple rehydration instructions
- Guaranteed 2nd generation and fully traceable to the original NCTC[®] and NCPF[®] source strain
- Certificates of Analysis available on our website
- Identification and Characterisation testing completed in accordance with our UKAS accredited scope (UKAS Testing Laboratory No. 2496)

Product data on non-selective agar

Code	Description
ABXXUL-1	Accu-Bac Ultra Low <100 CFU
ABXXL-1	Accu-Bac Low 50-150 CFU
ABXXM-1	Accu-Bac Medium 150-500 CFU
ABXXH-1	Accu-Bac High >500 CFU

